

Catalog Number

68504-15P

HPLC Column, RP C4 Butyl, 5um, 4.6mm ID x 150mm Length, 100A.
Cogent Brand. 1 EA.
MOST POPULAR.



Additional Info

This Cogent RP C4™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*butyl*) carbon ligand to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L26 Classification.

With a 5 µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Used for very hydrophobic or neutral compounds in Reversed Phase, this Butyl Phase is excellent for short analyses of small molecules. Expect Less retention than with a C8 column.

This is one of our Most Popular C4 Columns.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$848.36

Printed On 2024/05/19