

## Catalog Number

68504-SP100

**Semi Prep HPLC Column, RP C4, 5um, 100A. 10mm ID x 100mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP C4™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*butyl*) carbon ligand to the highest purity silica available.

With a 5 µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Used for very hydrophobic or neutral compounds in Reversed Phase, this Butyl Phase is excellent for short analyses of small molecules. Expect Less retention than with a C8 column.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,815.46**

Printed On 2024/05/19