

## Catalog Number

68508-15P-4

**HPLC Column, RP C8 Octyl, 5um, 4.0mm ID x 150mm Length, 100A.**  
Cogent Brand. 1 EA.  
*CUSTOM ORDER.*



### Additional Info

This Cogent RP C8™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 8 (*Octyl*) carbon silane to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L7 Classification.

With a 5 µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Used for hydrophobic or neutral compounds in Reversed Phase, this Octyl Phase is excellent for long chain fatty acids. Expect less retention than with a C18 column.

*This column is often chosen for most starting points of new methods or for Substitution for USP L7 Methods.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$848.36**

Printed On 2024/05/20