

## Catalog Number

68508-SP250

**Semi Prep HPLC Column, RP C8 Octyl, 5um, 100A. 10mm ID x 250mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP C8™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 8 (*Octyl*) carbon silane to the highest purity silica available.

The 5µm Particle works well with higher flow rates and produces lower back pressures than smaller particles. The diameter and long length of this column provides increased loading capacity compared with analytical columns, which is why it is chosen for purification of hydrophobic or neutral compounds, making it ideal for Semi Preparative HPLC.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,674.84**

Printed On 2024/05/19