

Catalog Number

68518-25P-4

HPLC Column, RP C18, 5um, 4.0mm ID x 250mm Length, 100A.

Cogent Brand. 1 EA.

CUSTOM ORDER.



Additional Info

This Cogent RP C18™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 18 carbon silane to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L1 Classification.

With a 5 µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids. Expect longer retention than with a C8 column.

This column is often chosen for starting points of new methods or for Substitution for USP L1 Methods.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$921.70

Printed On 2024/05/20