

## Catalog Number

68522-10P-2

**HPLC Column, RP PFP, 5um, 2.1mm ID x 100mm Length, 100A.**  
Cogent Brand. 1 EA.



### Additional Info

This Cogent RP PFP™ or RP Pentafluorophenyl Stationary Phase has been a premium value used by HPLC laboratories for many years and routinely called "F5". It is double end capped (*fully*) and bonded with a four carbon link to a benzene ring with 5 Fluorines; the highest purity silica available is used making this column ideal as an Alternate in USP Methods that require the L43 Classification.

With a 5µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Often chosen for samples with halogenated or polar functional groups.

*Packed and shipped in Reversed Phase Solvents.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$848.36**

Printed On 2024/05/19