

## Catalog Number

68524-P30-250

**Prep HPLC Column, RP Phenyl, 5um, 100A. 30mm ID x 250mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Cogent RP Phenyl™ Stationary Phase has been a premium value, and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available.

The 5µm Particle works well with higher flow rates.

The extra large diameter and long length of this column provides high loading capacity which is great for purification of isomers (positional), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, many Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds that have Phenyl Rings, making it ideal for Preparative HPLC.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$11,612.43**

Printed On 2024/05/20