

Catalog Number

70000-P21-150

Prep HPLC Column, Diamond Hydride, 4um, 100A. 21.2mm ID x 150mm Length. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Diamond Hydride Preparative, HPLC Column is made with 4um, 100A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for Purification or "Prep" of Peptides, Metabolites, Polar Compounds and other Analytes with amine groups when a large volume of sample is desired to be collected.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$6,746.39

Printed On 2024/05/20