

## Catalog Number

70200-15P-2

**HPLC Column, Diamond Hydride 2.0, 2.2um, 2.1mm ID x 150mm**  
Length, 120A. Cogent TYPE-C Brand. 1 EA.



### Additional Info

This Cogent Diamond Hydride 2.0 HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this smaller particle, silica hydride based stationary phase provides increased Efficiency, excellent Selectivity and faster Equilibration. Considered "Near UHPLC", this column can be used on any HPLC Instrument and still deliver improved efficiency.

Chosen for analysis of polar compounds that are difficult to retain in Reversed Phase, including polar and non-polar peptides, amino acids and other biological compounds such as metabolites. When using in Aqueous Normal Phase, this column is considered a Mass Spec Friendly Column.

*This column can be used for LCMS when increased resolution is desired.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,503.14**

Printed On 2024/05/20