

## Catalog Number

70200-HG2

**Guard Column Starter Set, Diamond Hydride 2.o**, Complete with 5 each 2.0mm ID x 10mm Long Cartridges with 2.2um, 120A and 1 each Guard Column Holder. Cogent TYPE-C Brand. 1/PK.



### Additional Info

Everything you need to get started using Cogent 2.o TYPE-C™ Guard Columns with your methods. This kit includes one "Universal" Stainless Steel [Guard Cartridge Column Holder](#) that can be used with any Cogent TYPE-C™ Guard Cartridge.

Five Cartridges, packed with Cogent Diamond Hydride 2.o Stationary Phase are included which will normally be sufficient protection for the lifetime of an analytical column, based on average use. Each Cartridge in this Kit is shipped in individual, recloseable Delrin Cases to keep them from drying out or being damaged.

*Use this Starter Kit with 2.1mm ID analytical columns packed with Diamond Hydride 2.o.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,478.80**

Printed On 2024/05/31