

Catalog Number

75001-P150

Prep HPLC Column, HPS Silica, 5um, 120A. 21.2mm ID x 150mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent HPS Silica™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. This proprietary Silica is not modified in any way and is ideal for purification of polar compounds by preparative HPLC.

With a 5 µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles and is often used for hydrophilic or polar compounds in Normal Phase HPLC.

Compare its performance to Inertsil and Luna products.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$9,470.76

Printed On 2024/05/20