

Catalog Number

75020-10P-2

HPLC Column, HPS Phenyl, 5um, 2.1mm ID x 100mm Length, 120A.
Cogent Brand. 1 EA.
MOST POPULAR.



Additional Info

This Cogent HPS Phenyl™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is fully end capped and is manufactured with a mono functional, 4 carbon ligand that has an attached *Phenyl* group and it is bonded to the highest purity, 120A silica available.

With a 5µm Particle, this Column works well with higher flow rates and lower back pressures compared with columns with smaller particles. It is often used to separate Isomers (*positional*), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds in Reversed Phase HPLC.

This column is often used with LCMS, compare its performance to Luna™ or Inertsil™ products.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,010.70

Printed On 2024/05/19