

## Catalog Number

75040-05P

**HPLC Column, HPS Amino, 5 $\mu$ m, 4.6mm ID x 50mm Length, 120A.**  
Cogent Brand. 1 EA.



### Additional Info

This Cogent HPS Amino™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using tri-functional bonding of a 3 carbon silane with a terminal amino group to our unique High Purity Silica. This Column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods.

With a 5  $\mu$ m Particle, this Column works well with higher flow rates and lower back pressures compared with smaller particles. Used for polar compounds with a negative charge such as Sugars and compounds with Carboxylic Acid Groups in Normal Phase or HILIC modes.

*Compare to Inertsil™ Amino and Luna™ Amino Columns. This column is sometimes used for scouting.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$857.64**

Printed On 2024/05/19