

Catalog Number

75040-15P

HPLC Column, HPS Amino, 5um, 4.6mm ID x 150mm Length, 120A.
Cogent Brand. 1 EA.
MOST POPULAR.



Additional Info

This Cogent HPS Amino™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Not end capped, it is manufactured using tri-functional bonding of a 3 carbon silane with a terminal amino group to our unique High Purity Silica. This Column is ideal as an Alternate in any USP Method that uses the L8 Classification or for new methods.

With a 5 µm Particle, this Column works well with higher flow rates and lower back pressures compared with smaller particles. Used for polar compounds with a negative charge such as Sugars and compounds with Carboxylic Acid Groups in Normal Phase or HILIC modes. Compare to Inertsil™ Amino and Luna™ Amino Columns.

This is one of our Most Popular HPS Amino™ Columns.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,060.16

Printed On 2024/05/19