

## Catalog Number

82560-P30-150

**Prep HPLC Column, NP Diol, 5um, 100A. 30mm ID x 150mm Length.**  
Cogent Brand. 1 EA.



### Additional Info

This Cogent NP Diol™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. Non end-capped, it is manufactured with the highest purity Silica available and bonded with a short hydrocarbon ligand that has a terminal, Diol Functional group.

The 5µm Particle in this Column works well with higher flow rates and yet produces acceptable back pressures. The large diameter and length of this column provides high loading capacity which is great for purification of compounds making it ideal for Preparative HPLC.

*Use for hydrophilic or polar compounds in HILIC or Normal Phase Prep HPLC instead of unmodified Silica Columns.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$7,917.08**

Printed On 2024/05/20