

## Catalog Number

84518-05P

**HPLC Column, RP C18 PE, 5um, 4.6mm ID x 50mm Length, 100A.**  
Cogent RP Brand. 1 EA.



### Additional Info

This Cogent RP C18 PE™ Stationary Phase has been a premium value and used by HPLC laboratories for many years.

It is an ideal Alternate for your USP Methods that use the L1 Classification and 100% Water as the Mobile Phase for hydrophobic samples. Very rugged, and often called "Polar Embedded", it can be used in typical Reversed Phase or *HILIC* modes.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids.

*This column is often chosen as a scouting column for preparative methods.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$743.63**

Printed On 2024/05/20