

Catalog Number

95001-0AP

Vial & Cap Kit includes 250 EA 1ml, Push Top, Amber, 8x40mm, Autosampler Vials and 250 EA matching Push Top Caps. AQ Brand. 250/PK.



Additional Info

1st hydrolytic, Type N51A Borosilicate Glass, amberized and considered low actinic. This product is manufactured in an ISO 9001:2015 Certified facility.

Push caps are made with ultra high purity polyethylene and have a thinned penetration point in a star configuration in the cap center to allow venting during injection and prevents the cap from grabbing the needle when it is extracted from the vials.

Use with Waters® WISP® 96 Position / Carousel autosamplers.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$70.95

Printed On 2024/05/19