

Catalog Number

95010-E-0AV-K

Vials, Headspace, Screw Top, Glass. Amber, 10ml, 18mm precision threads, 22.5x46mm outer dimensions, Rounded Bottom. CTC Compliant. BASIK Brand. (10 packs of 100) 1000/CS.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized and considered low actinic. Features 18mm precision threads for proper sealing.

Compatible with CTC®, LEAP Combi Pal®, Agilent® and other Headspace autosamplers. Not suitable for older versions of TurboMatrix 110® that were manufactured before 2006.

Lowest cost in class. Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$887.40

Printed On 2024/05/19