

## Catalog Number

95020-45S-W-K

**Caps, Headspace**, Screw Top, with White Silicone Rubber / Red PTFE Septa, fitted in Steel, Magnetic, Silver Caps with 8mm center holes. For use with 18mm Precision Thread Headspace Vials. MicroSolv Brand.



### Additional Info

Silicone Rubber / Ultraclean PTFE Septa has a durometer of 45° Shore A (hardness) and is 1.3mm thick. Cap has ND18 Screw Threads. Recommended with Agilent®, Perkin Elmer TurboMax® 16, 40 & 110 and other autosamplers.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,542.85**

Printed On 2024/05/19