


Catalog Number

95025-02-1S

Caps, Headspace, Crimp Top, with Grey Butyl Rubber Laminated PTFE Septa, fitted in Aluminum, Silver Caps with 10mm center holes. For use with 20mm Crimp Top Headspace Vials. MicroSolv Brand. 100/PK.



Additional Info

Butyl Rubber Septa that is laminated with PTFE on both sides and is 55° Shore A (hardness) and is 3.0mm thick. Temperature resistant from -40C to 120C.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$83.43

Printed On 2024/05/20