


## Catalog Number

95025-08-1S

**Caps, Headspace**, crimp top, with white silicone tubber / natural translucent PTFE, high purity septa, fitted in aluminum silver caps with 10 mm center holes. For use with 20 mm crimp top headspace vials.



### Additional Info

These high purity silicone rubber / PTFE septa are 3.16 mm (+/- 0.01 mm) thick, with a durometer of 45 Shore A, (hardness) and have 3.0 mil (0.0762 mm) thick "cast" PTFE liners. The ultra low bleed rubber that is made with a proprietary process for ultra durability, purity and low extractables.

Excellent choice for highly sensitive analyses for purity. The lower than standard durometer means that these are softer septa and therefore easier for penetrations by autosampler needles and have a

wide useful temperature range of -60° to 200°C.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$65.89**

Printed On 2024/05/20