


## Catalog Number

95025-14-3D

**Caps, Headspace**, Crimp Top, with Natural Silicone Rubber / Foil (125mil), High Purity Septa, fitted in Steel, Magnetic, Gold Caps with 5mm center holes. For use with 20mm Crimp Top Headspace Vials.



### Additional Info

Silicone Rubber / Aluminum Foil (125mil thick) Silver / Septa with Durometer of 50 Shore A, (hardness) and 3.0mm thick is good for high temperature analyses. No color has been added to this Septa. Often used with Perkin Elmer Instruments.

Wide temperature range from -60° to 220°C.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$169.88**

Printed On 2024/05/20