

## Catalog Number

95025-AG-06C

**Vials, Headspace**, Crimp Top, Glass. Clear, 6ml, 20mm Serum Finish, 22x38mm outer dimensions, Flat Bottom with a Long Neck. Agilent Compatible. BASIK Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Borosilicate Glass, Use with aluminum or steel crimp caps for proper sealing. Compatible with Agilent®, Carlo Erba®, Dani® and other Headspace autosamplers.

These vials are also used with Karl Fisher titrators by Mettler and Metrohm.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$63.28**

Printed On 2024/05/20