

## Catalog Number

95025-CRIMP-F

**Crimping Tool**, for 20mm, Pressure Relief, Aluminum Headspace Caps.  
Black handles. Hex adjustment key included. 1 EA.



### Additional Info

Chemical resistant coating makes these manual crimpers a great choice for laboratories.  
Adjustable height crimping jaws are made of a special, hardened alloy that produces a long life.

The height of the crimp is easily adjusted by screwing up or down the pressure block in the crimping head with the supplied hex key.

The maximum crimping pressure is controlled by the adjustable screw in the handle.

*Recommended for aluminum pressure relief (tear off seals) caps but not all steel crimp caps.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$503.33**

Printed On 2024/05/19