

Catalog Number

95025-E-AG20A

Vials, Headspace, Crimp Top, Glass. Amber, 20ml, 20mm Serum Finish, 22.5x75.5mm outer dimensions, Flat Bottom with a Long Neck. Agilent Compatible. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized and considered low actinic. Use with aluminum or steel crimp caps for proper sealing.

Compatible with Agilent® and other Headspace autosamplers

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$80.34

Printed On 2024/05/20