

## Catalog Number

9502C-0-0CV

**Vials, Crimp Top, Glass.** Clear, 2ml. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. MicroSolv Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass.

Use with 11mm aluminum crimp caps. These wide opening vials have *serum finish tops* that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$41.93**

Printed On 2024/05/19