


## Catalog Number

9502C-0E-0AV

**Vials, Crimp Top, Glass. Amber, 2ml. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. BASIK Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized. Amber vials are often used for light sensitive compounds and considered low actinic.

Use with 11mm aluminum crimp caps. These wide opening vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$33.58**

Printed On 2024/05/19