

## Catalog Number

9502C-0E-0DCV

**Vials, Crimp Top, Glass. Clear, 2ml. Deactivated. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Deactivated / Silanized Borosilicate Glass.

Use with 11mm aluminum crimp caps. These vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$46.20**

Printed On 2024/05/20