


Catalog Number

9502C-0E-WDAV

Vials, Crimp Top, Glass. Amber, 2ml, with Write-On Patch and fill lines. Deactivated. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Deactivated / Silanized Borosilicate Glass, amberized.

Amber vials are often used for light sensitive compounds and considered low actinic. Write-On Patch and fill lines are epoxy paint for reduced smearing with solvents.

Use with 11mm aluminum crimp caps. These wide opening vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$43.83

Printed On 2024/05/30