

Catalog Number

9502C-10CN

Caps, Snap Top, with Silicone Rubber / PTFE Septa, fitted in smooth, polyethylene, Natural Translucent Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

This cap is our original Snap Cap and has been on the market for 30 years. Click [HERE](#) for case quantities.

Natural Color Caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$40.33

Printed On 2024/05/20