

Catalog Number

9502C-10CN-M

Caps, Snap Top, with Silicone Rubber / PTFE Septa, fitted in smooth, polyethylene, Natural Translucent Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 1000/CS.



Additional Info

This cap is our original Snap Cap and has been on the market for 30 years. Click [HERE](#) for smaller package quantities..

Natural color caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$315.98

Printed On 2024/05/20