


## Catalog Number

9502C-30CG-M

**Caps, Snap Top**, with Silicone Rubber / PTFE Pre-Slit Septa, fitted in smooth, polyethylene, Green Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 1000/CS.



### Additional Info

This is one of the most popular snap caps on the market and has been available for over 25 years. *Pre-Slit septa are normally recommended for Waters Alliance instruments.*

Click [HERE](#) for smaller package quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$380.26**

Printed On 2024/06/01