


Catalog Number

9502C-34CT

Caps, Snap Top, with Silicone Rubber / PTFE Pre-Slit (Y Shape) Septa, fitted in smooth, polyethylene, Natural Transparent Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

Pre-Slit septa are normally recommended for Waters Alliance instruments.

Y Shape Pre-Slit offers a great target area for injector needles and are easier to penetrate but do not reseal as well as single slit septa.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$147.88

Printed On 2024/05/20