


Catalog Number

9502C-40P-S-M

Caps, Crimp Top, with PTFE Septa, fitted in aluminum, Sliver Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



Additional Info

This PTFE (only) Septa is for single injections or non-volatile samples as it will not reseal.

Buying in case quantities can save money and time. Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$233.80

Printed On 2024/05/20