

## Catalog Number

9502C-40T-S

**Caps, Crimp Top**, with Red Rubber / PTFE Septa, fitted in aluminum, Silver Caps. For use with 11mm Crimp or Snap Top, 12x32mm autosampler vials. BASIK Brand. 100/PK.



### Additional Info

Our most popular Crimp Cap. The Rubber is softer than Natural Red Rubber and provides better sealing after injection. Great for GC.

This is our most economical aluminum crimp cap, click [HERE](#) for case quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$17.48**

Printed On 2024/05/20