


Catalog Number

9502C-50CY-M

Caps, Snap Top, with White PTFE Septa, fitted in smooth, polyethylene, Yellow Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. MicroSolv Brand. 1000/CS.



Additional Info

This PTFE (only) Septa is for single injections or non-volatile samples as it will not reseal.
Buying in case quantities can save money and time.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$251.66

Printed On 2024/05/20