

Catalog Number

9502C-E-0DCV

Vials, Snap Top, Glass. Clear, 2ml. Deactivated. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Deactivated / Silanized Borosilicate Glass.

Lowest cost in class for budget minded laboratories.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$45.58

Printed On 2024/05/19