

Catalog Number

9502C-E-15C

Caps, Snap Top, with 0.04" thick, White Silicone Rubber / Natural PTFE Septa, fitted in smooth, polyethylene, Natural Transparent Caps. For use with 11mm, 12x32mm Snap Top autosampler vials. BASIK Brand.



Additional Info

Natural PTFE is very pure and does not have dyes that could extract from it.

Our most economical Silicone Rubber / PTFE Snap Cap. Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$28.59

Printed On 2024/05/20