

Catalog Number

9502C-PP

Vials, Snap Top, Plastic. Translucent, 1ml. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK.



Additional Info

These injection molded, translucent, clear plastic vials are made with virgin polypropylene and chosen for bio-analytical reasons or when glass vials are not indicated.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$52.60

Printed On 2024/05/19