

## Catalog Number

9502S-06-LCC

**Caps, Screw Top**, with Clear Silicone Rubber / PTFE, ultra-pure, LC-MS compatible LC Clean Septa, bonded in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials.



### Additional Info

Not Pre-Slit for better resealing properties. LC Clean™, LCMS compatible.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$52.73**

Printed On 2024/05/19