


Catalog Number

9502S-09-SP

Caps, Screw Top, Solid, with an ultra-clean, Silicone Rubber / PTFE liner, in "knurled", closed Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand 100/PK.



Additional Info

Solid top turns 9mm screw thread autosampler vials into storage vials. The liner has a silicone rubber backing to the sample facing PTFE which allows for a much better seal for volatile compounds or for solvents than PTFE only liners.

These caps are matched to work with our MicroSolv, AQ and RSA brand vials.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$58.20

Printed On 2024/05/19