

## Catalog Number

9502S-0AV-K

**Vials, Screw Top, Glass. Amber, 2ml, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. BASIK Brand. (10 packs of 100/PK). 1,000/CS.**



### Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass, amberized.

Features a wide mouth opening for easy filling and ND9 threads. Lowest cost in class.

Case quantity is the most economical value.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$310.63**

Printed On 2024/05/20