

Catalog Number

9502S-10-TEF

Caps, Screw Top, with White Silicone Rubber / Natural TEF Septa, fitted in smooth sided, Robo Style, soft polypropylene, Blue Caps. For use with 9-425 thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

White silicone rubber and white TEF fitted septa for a lower cost cap.

TEF is equally inert as PTFE for standard chromatography.

Compares to Thermo™ / National DP™ caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$36.30

Printed On 2024/05/19