

Catalog Number

9502S-11-MB-K

Caps, Screw Top, with Tan Silicone Rubber / White PTFE Septa, bonded in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQ Brand. 1000/CS.



Additional Info

Ultra-pure Septa and is a good choice for LCMS, GCMS & LCMS/MS. This cap is also recommended for Waters® Acquity™ autosamplers.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$416.41

Printed On 2024/05/19