

Catalog Number

9502S-31-MB-K

Caps, Screw Top, with Tan Silicone Rubber / PTFE Pre-Slit Septa, bonded in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials, AQ Brand. 1000/CS.



Additional Info

Pre-slit for easy penetration by autosampler needles and to prevent "Septa Plunge". Recommended for Waters® Alliance & Acquity instruments. Ultra-pure Silicone Rubber is a natural Tan color and will not extract pigments. Use this Cap for excellent results in LCMS/MS. Compare to "HPLC Certified" brands.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$629.54

Printed On 2024/05/19