


Catalog Number

9502S-32Y-M

Caps, Screw Top, with Silicone Rubber /PTFE Pre-Slit Septa, fitted in a smooth sided, Robo Style Yellow Caps. For use with 9-425 thread, 12x32mm autosampler vials. MicroSolv Brand. 1000/CS.



Additional Info

Pre-slit for easy penetration by autosampler needles and to prevent "Septa Plunge". Compares to Thermo / National Scientific DP® Cap.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$428.24

Printed On 2024/05/19