

Catalog Number

9502S-34C-B-M

Caps, Screw Top, with Silicone Rubber / PTFE **Partial Pre-Slit (Y Shape) Septa, fitted in "knurled" polypropylene, Blue Caps. For use with 9-425 thread, 12x32mm autosampler vials. Excellent for GC and GCMS.**



Additional Info

Partial slit, only in the rubber leaves PTFE to form a vapor barrier making this a good choice for volatile samples when an autosampler requires a pre-slit Septa.

The silicone rubber portion of the septa is 45 Shore and 1 mm (0.040 inch) thick. The Y shaped partial slit allows for better "resealing" of the rubber after initial injection.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$624.00

Printed On 2024/05/19