

Catalog Number

9502S-50C-N

Caps, Screw Top, with PTFE Septa, fitted in "knurled", polypropylene, Translucent Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



Additional Info

Rubber-free Cap is useful when rubber cannot be used with your analysis.

This PTFE (only) Septa is for single injections or non-volatile samples as it will not reseal

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$45.84

Printed On 2024/05/19