

Catalog Number

9502S-60C-B

Caps, Screw Top, with White Silicone Rubber / Clear Tefzel Septa, bonded in "knurled", polypropylene, Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK



Additional Info

Only 5 left in stock Inquire for discounts when purchased with AQ Brand Vials.

Need more? Click [HERE](#) for a replacement item.

Tefzel is inert to most solvents and this Cap is an economical alternative to Silicone Rubber / PTFE.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$33.91

Printed On 2024/05/19