

## Catalog Number

9502S-63C-B

**Caps, Screw Top**, with white silicone rubber / clear Tefzel pre-slit septa, bonded "knurled" polypropylene, blue caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK



### Additional Info

*Only 4 left in stock* Inquire for discounts when purchased with AQ Brand Vials.

Need more? Click [HERE](#) for a replacement item.

Tefzel is inert to most solvents and this Cap is an economical alternative to Silicone Rubber / PTFE. Pre-slit for easy penetration by autosampler needles. Recommended for Waters® Alliance.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$35.41**

Printed On 2024/05/19